

Product Name: CytoSep™ Cytology Funnel

Catalogue No.: M965-10FW

This document replaces any previous version

1. Product Description:

- **Cytology Funnel:** Disposable single funnel pre-assembled with white paper filter paper and cap included

2. Packaging:

- **Case:** 12 inner boxes of 40 / case of 480

3. Product Specifications:

- Compatible with the Shandon Cytospin® Cytocentrifuge:
- Material:
 - Funnel: Polypropylene
 - Cap: Polyethylene
 - Filter: White Cellulose filter board

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Paper conforms to 21CFR 176.170 & 21CFR 176.180
- Plastic conforms to FDA 21 CFR 177.1520
- **CE:** Product is CE marked
- **LATEX:** Material is Latex Free
- **BSE / TSE:** Material is BSE / TSE Free

5. Quality Assurance:

- No presence of contamination in plastic or paper
- Assembly and fit
- Visual attributes
- Dimension measurements

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store in at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination

8. Recommended Use:

- For use in **Shandon Cytospin® Cytocentrifuge** automated equipment, follow the equipment manufacturer's instructions for use

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*