

A 100

Upper basket

For the optimum loading of injector modules or inserts.



- Upper basket with two module connections
- for holding up to 2 injector modules or inserts
- Automatic closing of the connections when not in use
- Loading height with module with plastic supports 7.5 inches
- Clearances: width 19 1/4" (489 mm), length 17 5/16" (440 mm)

EAN: 4002515460610 / Material number: 09862320 / Old Material Number: 69610000D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
<b>Product category</b>	
Basket with module connection	•
Upper basket	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
<b>Application area</b>	
Reprocessing of pipettes/vials	•
Reprocessing of laboratory glassware	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	141
External dimensions, net width in inches	525
External dimensions, net depth in inches	528
External dimensions, gross height in inches	6 3/4 (170)
External dimensions, gross width in inches	21 7/8 (555)
External dimensions, gross depth in inches	22 1/4 (565)
Net weight in lbs	7 (2.80)
Gross weight in lbs	8 (3.49)