

TECHNICAL DATA SHEET

Part Number	934275B
Description	Tip 0,2 – 5 ml



Benefits :

- ✓ Eppendorf Biopur pipette tips in disposable Racks offer the highest biological purity grade
- ✓ Ensured PCR clean, sterile, pyrogen-, ATP- and DNA-free
- ✓ Sterility packaging is compliant to standards DIN EN ISO 11607 and DIN EN 868-2-10
- ✓ Non filtered
- ✓ Color : natural / rack violet
- ✓ Volume : 0,2 – 5 ml
- ✓ Length : 175 mm
- ✓ Packaging : 5 racks of 24 tips

Sales unit :

120 Tips