

## TECHNICAL DATA SHEET

---

<b>Part Number</b>	934479
<b>Description</b>	Tip 0,2 – 5 ml



---

### Benefits :

- ✓ EpT.I.P.S. Racks Eppendorf
- ✓ Long, narrow tips reach the base of narrow and high vessels without touching the top of the vessel edge
- ✓ Modern slim rack design with reclosable lid
- ✓ Non filtered
- ✓ Non sterile
- ✓ Color : natural / rack violet
- ✓ Volume : 0,2 – 5 ml
- ✓ Length : 175 mm
- ✓ Packaging : 5 racks of 24 tips

### Sales unit :

120 Tips

---