

MICRONIC | SAMPLE STORAGE | RACKS

MODULAR VERTICAL FREEZER RACK

DO YOU WANT A COST EFFECTIVE WAY TO STORE ANY STORAGE RACKS INTO A CHEST FREEZER? IF SO, THEN MICRONIC OFFERS YOU A FLEXIBLE AND LOW COST SOLUTION. THE MODULAR SYSTEM FITS IN ALL STANDARD STORAGE FREEZERS, IN TURN PROVIDING AN EFFECTIVE METHOD FOR STORING YOUR VALUABLE SAMPLES.

The Micronic Modular Vertical Freezer Rack provides a storage solution for samples into a chest freezer. Simply 'click' the individual extendable storage holders together to a tower of max. 15 holders.

One holder has a capacity of 1x Comorack-96 or 1x Roborack-96, using the 1.40ml tubes. When using the space saving 0.50ml or 0.75ml tubes in a Lobarack-96, the holder fits 2x Lobaracks-96, which doubles your storage capacity.

Using the Modular Vertical Freezer Rack system saves valuable freezer space.

For the secure tracking of your samples, all the tubes are available with a laser etched, 2D Data-Matrix or TraXis code on the bottom of the tube.

The 2D Data-Matrix or TraXis coded tubes can be decoded with one of our 2D code readers:

The Tracxer Code Reader can decode 96 tubes in 5 seconds, and the single tube scanner to simply decode one tube at a time.

For more information:
e-mail to sales@micronic.com

Micronic Europe B.V.
(Asia, Africa, Europe and Oceania)
Phone: +31 (0)320 277 070
Micronic America LLC
(North, Central and South America)
Phone: +1 484 480 3372

Please visit our website: www.Micronic.com



INFORMATION


- Cost effective solution
- Lightweight storage system
- Compatible with standard storage racks
- Easy to use modular system
- Class 7 clean-room production
- RNase/DNase and Pyrogene free




MICRONIC | SAMPLE STORAGE | RACKS

SPECIFICATION

| | |
|--------------|---|
| MATERIAL | ABS-PC |
| RESISTANCE | RESISTS HIGH CHEMICAL SOLVENTS |
| AUTOCLAVABLE | YES |
| CAPACITY | 1X COMORACK-96 OR ROBORACK-96 PER HOLDER 2X LOBORACK-96 PER HOLDER |



TECHNICAL DRAWING



ORDERING INFORMATION

| | | |
|--|-------------------------------|--------|
| ■ Modular Vertical Freezer Rack | | |
| MP34000 | Modular Vertical Freezer Rack | box/10 |



Micronic has a broader product range than the products listed under ordering information. For our complete product portfolio, please send your e-mail to sales@micronic.com or contact us by phone.