

MICRONIC | SAMPLE STORAGE | RACKS

LOBORACK-96

DO YOU REQUIRE A STORAGE RACK THAT CAN HOLD YOUR SPACE SAVING STORAGE TUBES? IF SO, THEN MICRONIC OFFERS THE LOW PROFILE LOBORACK-96 THAT CAN HOLD 96 INDIVIDUAL 0.50ML AND 0.75ML STORAGE TUBES.

The Loborack-96 can hold 96 individual tubes (8x12) in a 96 well configuration; individual retrieval of your highly valuable samples is easier than ever. This storage rack meets the global recognised standard footprint from the Society for Laboratory Automations and Screening (ANSI/SLAS), ensuring compatibility with most sample handling, storage systems and microplate storage racks.

The specially designed Loborack-96 fits the needs of automated systems:

- Space saving low profile design
- Unique Sure shot design ensures easy lead-in of tubes
- High chemical resistance
- Modular stackable lid and rack

The Loborack-96 can be closed securely with a transparent lid. With the unique and easy-to-use locking system, the lid stays perfectly closed throughout (long-term) storage and transportation.

The white Loborack-96 has a laser etched barcode that can never wear or fall off. The guaranteed unique 1D barcode ensures absolute sample traceability.


Available features:

- Sterile, gamma radiated
- Laser etched 1D barcode

For more information:
e-mail to sales@micronic.com

Micronic Europe B.V.
(Asia, Africa, Europe and Oceania)
Phone: +31 (0)320 277 070
Micronic America LLC
(North, Central and South America)
Phone: +1 484 480 3372

Please visit our website: www.Micronic.com



INFORMATION

- Industry standard Society for Laboratory Automations and Screening (ANSI/SLAS) footprint
- Secure and save locking system
- Absolute consistency
- Working temperature ultra low
- Available with custom specific barcodes
- Automation friendly
- Resists high chemical solvents
- Class 7 clean-room production
- RNase/DNase and Pyrogene free




MICRONIC | SAMPLE STORAGE | RACKS

SPECIFICATIONS

MATERIAL:	HIGHEST PURITY (MEDICAL GRADE) POLYPROPYLENE
ANSI/SLAS:	8X12 MICROPLATE STANDARD (9MM PITCH)
RESISTANCE:	RESISTS HIGH CHEMICAL SOLVENTS
AUTOCLAVABLE:	YES



TECHNICAL DRAWING



ORDERING INFORMATION

■ Lobarack-96		
MPW51004BC3	Lobarack-96 with cover (for 0.50ml tubes with Screw Cap or 0.75ml tubes with TPE or EVA Push Cap)	box/10
MPW51013BC3	Lobarack-96, with low cover (for 0.50ml tubes with TPE Caps)	box/10
MPW51015BC3	Lobarack-96, with high cover barcoded (for 0.75ml tubes with Screw Cap)	box/10



“Micronic Tubeselector”

Micronic has a broader product range than the products listed under ordering information. For our complete product portfolio, please send your e-mail to sales@micronic.com or contact us by phone.