

MICRONIC | SAMPLE STORAGE | CAPS

LOW PROFILE SCREW CAPS

DO YOU REQUIRE A SECURE SCREW CAP CLOSURE FOR YOUR SAMPLE STORAGE TUBES THAT SAVES ENORMOUS FREEZER SPACE DURING LONG-TERM SAMPLE STORAGE? IF SO, THEN MICRONIC OFFERS YOU THE LOW PROFILE SCREW CAP OPTIMIZED FOR THE MICRONIC SCREW CAP TUBES FROM 0.50ML TO 2.00ML IN A 96-WELL CONFIGURATION.

Micronic offers Low Profile Screw Caps that are compatible with the entire range of Micronic Screw Cap tubes in 96-well format: 0.50ml, 0.75ml, 1.10ml, 1.40ml and 2.00ml.

The Micronic Low Profile Screw Caps are ideal for biobanks and other laboratory facilities that need a space saving solution for their long-term sample storage at ultra-low temperatures. Since the tube height with Low Profile Screw Cap is now identical to the tube height with a TPE Cap (see picture below) - the same low profile rack cover can be used for both types of caps increasing operational flexibility for organizations using both types of storage tube closures.

The Low Profile Screw Caps are available in 12 different colors, ideal for your cap color coding practices (visual identification). The Silicon O-rings around the base of the Low Profile Screw Caps allow

for the tightest seal possible, and in turn provide maximum sample integrity for long-term storage.

The Micronic Screw Cap tubes are available non-coded, alphanumeric coded and 2D coded (Data-Matrix or TraXis), and manufactured from the highest-purity polypropylene. The 2D codes are encrypted with a laser in the polypropylene code area, so that they are guaranteed resistant to solvents and ensuring a reliable code reading.

Available features

- Sterile, gamma radiated
- Precapped on Micronic Screw Cap tubes

For more information:
e-mail to sales@micronic.com



Micronic Europe B.V.
(Asia, Africa, Europe and Oceania)
Phone: +31 (0)320 277 070

Micronic America LLC
(North, Central and South America)
Phone: +1 484 480 3372

Please visit our website: www.Micronic.com

INFORMATION



- Ideal for long-term sample storage at ultra-low temperatures
- Space saving solution
- Torx10 screw head
- Silicon O-ring for optimal closure
- Available in 12 different colors
- Automation friendly
- Clean-room Class 7 production
- RNase/DNase and Pyrogene free



MICRONIC | SAMPLE STORAGE | CAPS

SPECIFICATIONS

MATERIAL	PBTP WITH SILICON O-RING
RESISTANCE	CHEMICAL SOLVENTS
COLORS	GREY, WHITE, YELLOW, ORANGE, RED, PINK, PURPLE, BLUE, LIGHT BLUE, LIGHT GREEN, GREEN AND BLACK
AUTOCLAVABLE TEMPERATURE	UP TO TEMPERATURES OF +121°C VAPOR PHASE LN2 TO +121°C



ORDERING INFORMATION

■ Screw Caps Low Profile for tubes in 96-well format

MP53320	Screw Caps Low Profile for tubes in 96-well format - Grey	box/960
MP53301	Screw Caps Low Profile for tubes in 96-well format - White	box/960
MP53304	Screw Caps Low Profile for tubes in 96-well format - Yellow	box/960
MP53300	Screw Caps Low Profile for tubes in 96-well format - Orange	box/960
MP53302	Screw Caps Low Profile for tubes in 96-well format - Red	box/960
MP53321	Screw Caps Low Profile for tubes in 96-well format - Pink	box/960
MP53309	Screw Caps Low Profile for tubes in 96-well format - Purple	box/960
MP53303	Screw Caps Low Profile for tubes in 96-well format - Blue	box/960
MP53306	Screw Caps Low Profile for tubes in 96-well format - Light Blue	box/960
MP53323	Screw Caps Low Profile for tubes in 96-well format - Light Green	box/960
MP53305	Screw Caps Low Profile for tubes in 96-well format - Green	box/960
MP53322	Screw Caps Low Profile for tubes in 96-well format - Black	box/960



Micronic has a broader product range than the products listed under ordering information. For our complete product portfolio, please send your e-mail to sales@micronic.com or contact us by phone.