



Bel-Art Products, a part of SP
661 Route 23 South, Wayne, New Jersey 07470
Tel: 973-694-0500 * Fax: 973-694-7199
www.sp-belart.com

CERTIFICATE OF ANALYSIS

This is to verify that all items noted below are in conformance with the Purchase Order Contract, drawings, specifications, and other applicable documentation; that all process certifications, material certifications, are on file at this or another Bel-Art facility and are available for review by the customer and or Government Representatives upon request.

Product #: 42045-0500 **Mfg. Date:** N/A **Lot #:** GENERAL

Description: Desiccant in a Bag, 500g

Material: Microfiber cloth, Indicating Silica Gel Beads with Cobalt Dichloride

Dimensions: 6" (15cm) L x 2" (5cm) W x 8.25" (21cm) H, 500g

Other: Beads turn pink from blue when saturated with moisture. Can be regenerated in a vented oven at 150°C (300°F) for 3 hours. Silica Gel is inert, non-corrosive, and non-toxic.
Not FDA Grade. Absorbs approximately 20% moisture by wt. Moisture content- max: 4.0%.
No Shelf Life if desiccant remains in sealed foil pouch, suggested to be used within 3 years.

Manufactured by/for: Bel-Art Products, a part of SP
2024 Broad Street Pocomoke City, MD 21851
An ISO-9001:2015 facility (certificate US19/819943409)
Made in China

Reviewed and Approved By:

Signature:

Printed Name: Nancy K. Trogani

Title: Technical Support Analyst

Date: 6/17/2022