


TECHNICAL DATA SHEET

Part Nbr	949001
Designation	Reusable rack



Description :

- ✓ Rack designed for Labcon 15 ml centrifuge tubes
- ✓ Autoclavable (122°C)
- ✓ Freezable (-80°C)
- ✓ Can be stacked together to minimize storage space
- ✓ Interlocking design so you can connect more than one rack and create a handy work station
- ✓ Color : white
- ✓ 1 rack is designed to hold 25 tubes of 15ml
- ✓ High of the rack 31 mm

Sales unit :

5 Reusable racks