

Product Description

Catalog Number: 3095-870-008

Product Description: 50ml Centrifuge Tube. Free-standing
Polypropylene. Molded graduations. Flat caps

Product Dimensions:

Length: 115mm

Inside Diameter: 25mm

Outside Diameter: 29mm

Packaging: Bulk - 500 tubes per bag. (Tubes and caps in separate bags)

Can be spun to 10,000 RCF.

MicroBiological testing:

Certified testing is performed lot by lot. These tests include Nuclease(Rnase Dnase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.