

## Product Description

**Catalog Number:** 1030-800-000

**Product Description:** 1-200ul yellow pipet tip. Polypropylene. Non-sterile.  
Loose pack.

**Product Dimensions:**

Length of Tip: 1.955" (49.657mm)  
Outside Diameter at tip: .030" (.762mm)  
Inside Diameter at tip: .020" (.508mm)  
Outside Diameter at opening: .30" (7.62mm)  
Inside Diameter at opening: .21" (5.33mm)

**Packaging:** 1000 tips per pack  
10 packs per case (10,000 tips total)

Autoclavable.

**MicroBiological testing:**

Certified testing is performed lot by lot. These tests include Nuclease(Rnase Dnase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

**Performance Testing:**

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.