

Product Description

Catalog Number: 1093-800-000-9

Product Description: 1-200 μ l pipet tip with Bevel Point
Polypropylene. Clear. Non-Sterile

Product Dimensions:

Length of Tip: 1.96" (49.78mm)
Outer Diameter at dispensing end: .03" (.76mm)
Inner Diameter at dispensing end: .02" (.50mm)
Outer Diameter at pipettor: .30" (7.62mm)
Inner Diameter at pipettor: .21" (5.33mm)

Packaging: Bulk (loose)
1000 tips per bag/pack
10 bags/packs per case (10,000 tips total)

Autoclavable.

MicroBiological testing:

Certified testing is performed lot by lot. These tests include Nuclease (RNAse and DNAse), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.