

Product Description

Catalog Number: 1045-800-000

Product Description: 100-1250ul pipet tip. Extended length with Flextop and UltraFine point. Polypropylene. Non-sterile. Loose pack.

Product Dimensions:

Length of Tip: 3.500" (mm)

Outside Diameter at tip: .036" (.914mm)

Inside Diameter at tip: .020" (.508mm)

Outside Diameter at opening: .290" (7.36mm)

Inside Diameter at opening: .242" (6.14mm)

Packaging: 1000 tips per pack.
10 packs per case (10,000 tips total)

Autoclavable.

MicroBiological testing:

Certified testing is performed lot by lot. These tests include Nuclease(Rnase Dnase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.