



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Eclipse™ FlexTop™ 1250 μ L Wide Orifice Pipet Tips, in Resealable Bags

Fluid Measure	μ L	Freezable	NA
Fluid Volume	1250	Freeze Temp	
Height (mm)	88.9	Autoclavable	121°C 15min
Width (mm)		Sterile	No
Length (mm)		Sterile Processing	NA
Pack Qty	1000	Centrifugation	NA
Case Qty	10000	Centrifugation Speed (RCF)	
Package Type	Resealable Bags	Shelf Life (Months)	48
Material	Polypropylene	Certification	
Cap Material		Regulatory	

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing: ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)