


# Technical Product Data

Technical Data Sheet

Labcon Product Number 3191-335-018-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

## 50 mL SuperClear® Centrifuge Tubes, 25 per Rack, Sterile

Fluid Measure	mL	Freezable	Yes
Fluid Volume	50	Freeze Temp	-80°C
Height (mm)	115	Autoclavable	121°C 15min
Width (mm)	29	Sterile	Yes
Length (mm)		Sterile Processing	Certified as sterile with a sterility assurance level (SAL) of 10 <sup>-6</sup>
Pack Qty	50	Centrifugation	Yes
Case Qty	500	Centrifugation Speed (RCF)	20000
Package Type	Cardboard Fiber Rack	Shelf Life (Months)	60
Material	Polypropylene	Certification	ISO 9001
Cap Material	Polypropylene	Regulatory	U.S. FDA 21 CFR

**Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.**

### Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

IATA 95kPa certified to conform with USDOT HMR, 49 C.F.R., Parts 171-180 and IATA DGR UN3373 PI650 mailing tube standards and packaging instructions for diagnostic and infectious specimens.