



Labcon, North America  
3700 Lakeville Highway  
Suite 200  
Petaluma, CA 94954-5671 U.S.A.  
707.766.2100  
Fax 707.766.2199  
www.labcon.com

## Product Description

**Catalog Number:** 3191-870-008-9

**Product Description:** 50mL SuperClear Centrifuge Tube  
Polypropylene. Non-Sterile. Flat cap.  
Printed Graduations. Anti-Smear Writing Panel

**Product Dimensions:** Length of Tube: 4.51" (114.55mm)  
Outside Diameter: 1.16" (29.46mm)  
Inside Diameter: 1.08" (27.43mm)

**Packaging:** Bulk – 500 tubes per case (500 caps in separate bags)

Autoclavable to 122 C. Freezable. Tubes can be centrifuged to 20,000 RCF.

**MicroBiological testing:**

Certified testing is performed lot by lot. These tests include Nuclease (RNase and DNase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

**Performance Testing:**

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.