



Labcon, North America
3700 Lakeville Highway
Suite 200
Petaluma, CA 94954-5671 U.S.A.
707.766.2100
Fax 707.766.2199
www.labcon.com

Product Description

Catalog Number: 3141-335-018-9

Product Description: 15mL SuperClear Centrifuge Tube
Polypropylene. Sterile. Graduations.
Polyethylene Flat Cap Cardboard Fiber Rack

Product Dimensions:

Length of Tube: 4.69" (119.12mm)
Outside Diameter: .69" (17.52mm)
Inside Diameter: .58" (14.83mm)

Packaging: Racked (with attached caps)
25 tubes per rack.
2 racks per pack (50 tubes total)
10 packs per case (500 tubes total)

Autoclavable. Freezable. Tubes can be centrifuged to 17,000 RCF.

MicroBiological testing:

Certified testing is performed lot by lot. These tests include Nuclease (RNAse and DNAse), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Sterility:

Certified as sterile with a sterility assurance level (SAL) of 10^{-6}