



Labcon, North America
3700 Lakeville Highway
Suite 200
Petaluma, CA 94954-5671 U.S.A.
707.766.2100
Fax 707.766.2199
www.labcon.com

Product Description

Catalog Number: 1059-165-000

Product Description: 1-200uL pipet tip. UltraFine Point Polypropylene. Filtered. Sterile.

Product Dimensions:

Length of Tip: 2.05" (52.07mm)

Outside Diameter at tip: .050" (1.016mm)

Inside Diameter at tip: .02" (.508mm)

Outside Diameter at opening: .288" (7.315mm)

Inside Diameter at opening: .21" (5.334mm)

Packaging: Individually-wrapped
200 tips per pack
10 packs per case (2000 tips total)

Autoclavable.

MicroBiological testing:

Certified testing is performed lot by lot. These tests include Nuclease(Rnase Dnase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

Sterility Testing: Those products designated as Sterile are irradiated by gamma radiation. The validated dosage level for these products is in a range between 12-25kGy.