



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

225 mL PerformR® Centrifuge Tubes, in Bags, Sterile

Fluid Measure	mL	Freezable	Yes
Fluid Volume	225	Freeze Temp	-80°C
Height (mm)	137	Autoclavable	121°C 15min
Width (mm)	61	Sterile	Yes
Length (mm)		Sterile Processing	Certified as sterile with a sterility assurance level (SAL) of 10 ⁻⁶
Pack Qty	2	Centrifugation	Yes
Case Qty	48	Centrifugation Speed (RCF)	8000
Package Type	Bag	Shelf Life (Months)	36
Material	Polypropylene	Certification	
Cap Material	Polypropylene	Regulatory	

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing: ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)