



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

15 mL ProtectR® Dry Ice Storage Tubes in IntegraPack®, 10 per Bag

Fluid Measure	mL	Freezable	Yes
Fluid Volume	15	Freeze Temp	-80°C
Height (mm)	118	Autoclavable	121°C 15min
Width (mm)	17	Sterile	No
Length (mm)		Sterile Processing	NA
Pack Qty	50	Centrifugation	Yes
Case Qty	100	Centrifugation Speed (RCF)	12500
Package Type	IntegraPack®	Shelf Life (Months)	48
Material	Polypropylene	Certification	
Cap Material	Polypropylene	Regulatory	

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing: ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

IATA 95kPa certified to conform with USDOT HMR, 49 C.F.R., Parts 171-180 and IATA DGR UN3373 PI650 mailing tube standards and packaging instructions for diagnostic and infectious specimens.

