

Dimensions are inches. Do not scale




Image is not actual size. Colors may differ.

## 50 mL ProtectR® Dry Ice Storage Tubes in IntegraPack®, 10 per Bag, Sterile

Fluid Measure	mL	Freezable	Yes
Fluid Volume	50	Freeze Temp	-80°C
Height (mm)	115	Autoclavable	121°C 15min
Width (mm)	29	Sterile	Yes
Length (mm)		Sterile Processing	Certified as sterile with a sterility assurance level (SAL) of 10 <sup>-6</sup>
Pack Qty	50	Centrifugation	Yes
Case Qty	100	Centrifugation Speed (RCF)	12500
Package Type	IntegraPack®	Shelf Life (Months)	36
Material	Polypropylene	Certification	
Cap Material	Polypropylene	Regulatory	

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

**Additional Testing:** ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

IATA 95kPa certified to conform with USDOT HMR, 49 C.F.R., Parts 171-180 and IATA DGR UN3373 PI650 mailing tube standards and packaging instructions for diagnostic and infectious specimens.

