



PerformR® Centrifuge Tubes

For your everyday centrifugation and storage needs



Reliable performance • Higher spin rate • High quality medical-grade resin

www.labcon.com



PerformR® Centrifuge Tubes are designed for your everyday centrifugation and storage needs. These tubes include large white solvent resistant labeling areas and black easy-to-read graduations. Available in packaging styles for all your needs, including reusable polypropylene racks, sterile-bulk bags, and non-sterile bulk bags. Always use appropriate cushion or adapters and biohazard containment precautions.

VISIBILITY PerformR® tubes are made from clear, medical-grade polypropylene resin that allows you see sample volume and color easily. These tubes feature printed graduations from UV cured inks that won't rub off.

RELIABILITY PerformR® tubes are developed using our own in-house laboratory and tested with everyday basic procedures in real laboratories. Printed graduations are tested for accuracy during every production lot.


SECURITY PerformR® tubes are available with either a flat cap that is easy to label or a plug style cap for horizontal applications. Both caps offer a secure seal.

SUSTAINABILITY PerformR® tubes are designed for maximum product recyclability and manufactured in a solar-powered facility with unique waste-free molding machines. Tube graphics and graduations are printed with UV cured inks which significantly reduces the emissions of harmful volatile organic compounds (VOCs). Offered in reusable and recyclable plastic racks.


PerformR® tubes are available non-sterile or sterile. Our sterile products are radiation sterilized following a validated ISO 11137 method.

Special Features

Flat Or Plug Style Caps



Printed Graduations




Solvent-Resistant Writing Panel



Recyclable



Spins To 12,500 RCF



Polypropylene tubes and can be autoclaved (122°C) or frozen (-80°C).

PerformR Centrifuge Tube Cap Styles:



FLAT CAPS

- Secure seal
- Easy to label
- High-density polyethylene
- Autoclavable* (122°C)
*with caps loosely attached to centrifuge tubes



PLUG CAPS

- Secure seal
- Ideal for extended storage
- High-density polyethylene
- Autoclavable* (122°C)
*with caps loosely attached to centrifuge tubes

PACKAGING OPTIONS



PLASTIC RACKS

25 Tubes Per Rack
Sterile

- Reusable and recyclable
- Stackable and interlockable to form benchtop workstations
- Numbered wells
- Autoclavable (122°C)
- Freezable (-80°C)



BAGS

50 Tubes Per Bag
Sterile

- Bags have a validated secure seal for purity
- Great value
- Economical choice



BULK

500 Tubes Per Box

Tubes and caps in separate bags • Excellent economical choice

ORDERING INFORMATION



15mL PerformR® Centrifuge Tubes

Product #	Cap	Sterile	Packaging	Quantity Per Case
3146-335-000-9	Flat	✓	Plastic Rack	12 Racks (300 tubes)
3136-345-000-9	Flat	✓	Bag	10 Bags (500 tubes)
3146-870-000-9	Flat	-	Bulk	1 Box (500 tubes)
3147-335-000-9	Plug	✓	Plastic Rack	12 Racks (300 tubes)
3137-345-000-9	Plug	✓	Bag	10 Bags (500 tubes)
3147-870-000-9	Plug	-	Bulk	1 Box (500 tubes)

Note: 15mL tubes don't work with Thermo (010-0378) centrifuge adapters



50mL PerformR® Centrifuge Tubes

Product #	Cap	Sterile	Packaging	Quantity Per Case
3196-335-000-9	Flat	✓	Plastic Rack	12 Racks (300 tubes)
3186-345-000-9	Flat	✓	Bag	10 Bags (500 tubes)
3196-870-000-9	Flat	-	Bulk	1 Box (500 tubes)
3197-335-000-9	Plug	✓	Plastic Rack	12 Racks (300 tubes)
3187-345-000-9	Plug	✓	Bag	10 Bags (500 tubes)
3197-870-000-9	Plug	-	Bulk	1 Box (500 tubes)

labcon®



MADE IN USA

+1.707.766.2100 | info@labcon.com | www.labcon.com

Quality & Compliance



Product details are current when published and subject to change without notice. Labcon®, the Labcon logo, and variations on the foregoing are trademarks, or registered trademarks of Labcon North America. © Copyright 2019 Labcon. All rights reserved.

REV 0_2.25.19