



Labcon, North America
3700 Lakeville Highway
Suite 200
Petaluma, CA 94954-5671 U.S.A.
707.766.2100
Fax 707.766.2199
www.labcon.com

Product Description

Catalog Number: 3131-345-008-9

Product Description: 15ml Centrifuge Tube with Flat Cap*
Polypropylene. Sterile Graduated.

Product Dimensions:

Length of Tube: 4.71" (119.6mm)
Outer Diameter: .68" (17.3mm)
Inner Diameter: .60" (15.2mm)

Packaging: Bagged
50 tubes per bag.
10 bags per case (500 tubes total)

Freezable. Tubes can be centrifuged to 17,000 RCF.

Autoclavable. (Loosen the caps of tubes to prevent collapse when autoclaving)

*Caps include an elastomeric molded-in sealing ring

MicroBiological testing:

Certified testing is performed lot by lot. These tests include Nuclease (RNase DNase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Sterility:

Certified as sterile with a sterility assurance level (SAL) of 10^{-6}