



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

15 mL PerformR® Polystyrene Centrifuge Tubes, in Bags, Sterile

Fluid Measure	mL	Freezable	No
Fluid Volume	15	Freeze Temp	-80°C
Height (mm)	118	Autoclavable	No
Width (mm)	17	Sterile	Yes
Length (mm)		Sterile Processing	Certified as sterile with a sterility assurance level (SAL) of 10 ⁻⁶
Pack Qty	50	Centrifugation	Yes
Case Qty	500	Centrifugation Speed (RCF)	3000
Package Type	Bag	Shelf Life (Months)	36
Material	Polystyrene	Certification	
Cap Material	Polypropylene/HDPE Liner	Regulatory	

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing: ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

DOMINIQUE DUTSCHER SAS