



Labcon, North America
3700 Lakeville Highway
Suite 200
Petaluma, CA 94954-5671 U.S.A.
707.766.2100
Fax 707.766.2199
www.labcon.com

Product Description

Catalog Number: 3131-325-008-9

Product Description: 15ml Centrifuge Tube with printed graduations
Polypropylene. Sterile. With attached PP Flat Caps.
ViewPoint writing panel

Product Dimensions:

Length of Tube: 4.69" (119.15mm)
Outside Diameter: .69" (17.52mm)
Inside Diameter: .58" (14.83mm)

Packaging: Integrapack
50 tubes per pack (5 sub-packs of 10 tubes)
10 packs per case (500 tubes total)

Autoclavable. Freezable.
Tubes can be centrifuged to 17,000 RCF.

MicroBiological testing:

Certified testing is performed lot by lot. These tests include Nuclease (RNase and DNase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Sterility:

Certified as sterile with a sterility assurance level (SAL) of 10^{-6}