



Labcon, North America  
3700 Lakeville Highway  
Suite 200  
Petaluma, CA 94954-5671 U.S.A.  
707.766.2100  
Fax 707.766.2199  
www.labcon.com

## Product Description

**Catalog Number:** 3132-340-008-9

**Product Description:** 15ml SuperClear Centrifuge Tube,  
Polypropylene. Non-Sterile. Plug Style cap.

**Product Dimensions:**

Length of Tube: 4.71" (119.6mm)  
Outside Diameter: .68" (17.3mm)  
Inside Diameter: .60" (15.2mm)

**Packaging:** 50 tubes with attached caps per bag.  
10 bags per case (500 tubes total)

Freezable. Autoclavable. Centrifuge to 17,000 RCF.

**MicroBiological testing:**

Certified testing is performed lot by lot. These tests include Nuclease (RNase and DNase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

**Performance Testing:**

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.