



Labcon, North America
3700 Lakeville Highway
Suite 200
Petaluma, CA 94954-5671 U.S.A.
707.766.2100
Fax 707.766.2199
www.labcon.com

Product Description

Catalog Number: 3132-345-008-9

Product Description: 15ml SuperClear Centrifuge Tube,
Polypropylene. Sterile. Plug Style cap.

Product Dimensions:

Length of Tube: 4.71" (119.6mm)

Outside Diameter: .68" (17.3mm)

Inside Diameter: .60" (15.2mm)

Packaging: 50 tubes with attached caps per bag.
10 bags per case (500 tubes total)

Freezable. Autoclavable. Centrifuge to 17,000 RCF.

MicroBiological testing:

Certified testing is performed lot by lot. These tests include Nuclease (RNase and DNase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Sterility:

Certified as sterile with a sterility assurance level (SAL) of 10^{-6} .