



Labcon, North America
3700 Lakeville Highway
Suite 200
Petaluma, CA 94954-5671 U.S.A.
707.766.2100
Fax 707.766.2199
www.labcon.com

Product Description

Catalog Number: 3134-345-001-9

Product Description: 15ml Metal Free Centrifuge Tube
Polypropylene. Sterile. Graduations
Flat natural-colored polyethylene cap

Product Dimensions:

Length of Tube: 4.69" (119.12mm)
Outer Diameter: .69" (17.52mm)
Inner Diameter: .58" (14.73mm)

Packaging: 50 tubes per bag.
10 bags per case (500 tubes total)

Autoclavable. Freezable. Tubes can be centrifuged to 17,000 RCF.

MicroBiological testing:

Certified testing is performed lot by lot. These tests include Nuclease (RNAse and DNAse), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.