

Product Description

Catalog Number: 3140-335-028-9

Product Description: 15mL Centrifuge Tube with etched graduations
Polystyrene. Sterile. Polyethylene Flat cap*

Product Dimensions:

Length of Tube: 4.69" (119.12mm)
Outside Diameter: .69" (17.52mm)
Inside Diameter: .58" (14.73mm)

Packaging: Racked (cardboard fiber rack)
50 tubes per pack (2 racks of 25 tubes each)
10 packs per case (500 tubes total)

Tubes can be centrifuged to 3,000 RCF.

* High Density Polyethylene cap with foamed polyethylene liner.

MicroBiological testing:

Certified testing is performed lot by lot. These tests include Nuclease (RNase and DNase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Sterility:

Certified as sterile with a sterility assurance level (SAL) of 10^{-6}