



Labcon, North America  
3700 Lakeville Highway  
Suite 200  
Petaluma, CA 94954-5671 U.S.A.  
707.766.2100  
Fax 707.766.2199  
www.labcon.com

## Product Description

**Catalog Number:** 3141-335-028-9

**Product Description:** 15mL SuperClear Centrifuge Tube  
Polypropylene. Sterile. Graduations.  
Polyethylene Flat Cap Plastic Rack

**Product Dimensions:**

Length of Tube: 4.69" (119.12mm)  
Outside Diameter: .69" (17.52mm)  
Inside Diameter: .58" (14.83mm)

**Packaging:** Racked (with attached caps)  
25 tubes per rack.  
2 racks per pack (50 tubes total)  
10 packs per case (500 tubes total)

Autoclavable. Freezable. Tubes can be centrifuged to 17,000 RCF.

**MicroBiological testing:**

Certified testing is performed lot by lot. These tests include Nuclease (RNase and DNase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

**Performance Testing:**

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

**Sterility:**

Certified as sterile with a sterility assurance level (SAL) of  $10^{-6}$