



Labcon, North America  
3700 Lakeville Highway  
Suite 200  
Petaluma, CA 94954-5671 U.S.A.  
707.766.2100  
Fax 707.766.2199  
www.labcon.com

## Product Description

**Catalog Number:** 3181-165-008-9

**Product Description:** 50mL SuperClear Centrifuge tube  
Printed graduations. Sterile. Polypropylene.  
Polyethylene flat cap.

**Product Dimensions:**  
Length of Tube: 4.51" (114.55mm)  
Outer Diameter of Tube 1.16" (29.46mm)  
Inner Diameter of Tube 1.07" (27.17 mm)

**Packaging:** Individually wrapped. 100 tubes per case.

**MicroBiological testing::**  
Certified testing is performed lot by lot. These tests include Nuclease (RNase and DNase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

**Performance Testing:**  
Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

**Sterility:**  
Certified as sterile with a sterility assurance level (SAL) of  $10^{-6}$