



Labcon, North America
3700 Lakeville Highway
Suite 200
Petaluma, CA 94954-5671 U.S.A.
707.766.2100
Fax 707.766.2199
www.labcon.com

Product Description

Catalog Number: 3181-325-008-9

Product Description: 50mL SuperClear Centrifuge tube
Printed graduations. Sterile. Polypropylene.
Polyethylene flat cap.

Product Dimensions:

Length of Tube: 4.51" (114.55mm)
Outer Diameter of Tube 1.16" (29.46mm)
Inner Diameter of Tube 1.07" (27.17 mm)

Packaging:

“IntegraPack”
5 packs of 10 tubes per pack
10 packs per case (500 tubes total)

MicroBiological testing::

Certified testing is performed lot by lot. These tests include Nuclease (RNase and DNase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Sterility:

Certified as sterile with a sterility assurance level (SAL) of 10^{-6}