



Labcon, North America
3700 Lakeville Highway
Suite 200
Petaluma, CA 94954-5671 U.S.A.
707.766.2100
Fax 707.766.2199
www.labcon.com

Product Description

Catalog Number: 3181-340-008-9

Product Description: 50mL SuperClear® Centrifuge Tubes with printed graduations, anti-smear writing panel and SuperSeal™ flat caps with integral elastomeric seal.

Tube Material: Polypropylene

Cap Material: Polypropylene

Product Dimensions:

Length of Tube: 4.51" (114.6mm)

Outer Diameter: 1.16" (29.5mm)

Inner Diameter: 1.08" (27.4mm)

Packaging: Bagged
50 tubes per bag
10 bags per case (500 tubes total)

Freezable. Tubes can be centrifuged to 20,000 RCF.

Autoclavable. (Remember to loosen the caps of tubes to prevent collapse when autoclaving)

MicroBiological testing:

Certified testing is performed lot by lot. These tests include Nuclease (RNase and DNase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.