



Labcon, North America  
3700 Lakeville Highway  
Suite 200  
Petaluma, CA 94954-5671 U.S.A.  
707.766.2100  
Fax 707.766.2199  
www.labcon.com

## Product Description

**Catalog Number:** 3191-335-028-9

**Product Description:** 50ml SuperClear Centrifuge Tube  
Polypropylene. Sterile. Graduations.  
Polyethylene Flat Cap. Plastic Rack

**Product Dimensions:**

Length of Tube: 4.51" (114.55mm)

Outer Diameter: 1.16" (29.46mm)

Inner Diameter: 1.08" (27.43mm)

**Packaging:** Racked (with attached caps)  
25 tubes per rack  
2 racks per pack (50 tubes total)  
10 racks per case (500 tubes total)

Autoclavable. Freezable. Tubes can be centrifuged to 20,000 RCF.

**MicroBiological testing:**

Certified testing is performed lot by lot. These tests include Nuclease (RNAse and DNAse), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

**Performance Testing:**

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

**Sterility:**

Certified as sterile with a sterility assurance level (SAL) of  $10^{-6}$