



Labcon, North America  
3700 Lakeville Highway  
Suite 200  
Petaluma, CA 94954-5671 U.S.A.  
707.766.2100  
Fax 707.766.2199  
www.labcon.com

## Product Description

**Catalog Number:** 3192-330-018-9

**Product Description:** 50mL SuperClear Centrifuge Tube  
Polypropylene Non-sterile. Printed Graduations  
Plug-Style Polyethylene Cap

**Product Dimensions:**

Length of Tube: 4.51" (114.55mm)  
Outer Diameter of Tube 1.16" (29.46mm)  
Inner Diameter of Tube 1.08" (27.43mm)

**Packaging:** Racked (cardboard fiber racks)  
25 tips per rack (with attached caps)  
2 racks per pack (50 tubes)  
10 packs per case (500 tubes total)  
500 tubes/caps

Autoclavable. Freezable. Tubes can be centrifuged to 20,000 RCF.

**MicroBiological testing:**

Certified testing is performed lot by lot. These tests include Nuclease (RNase and DNase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

**Performance Testing:**

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.