



Labcon, North America
3700 Lakeville Highway
Suite 200
Petaluma, CA 94954-5671 U.S.A.
707.766.2100
Fax 707.766.2199
www.labcon.com

Product Description

Catalog Number: 3192-335-018-9

Product Description: 50mL SuperClear Centrifuge Tube
Polypropylene Sterile. Printed Graduations
Plug-Style Polyethylene Cap

Product Dimensions:

Length of Tube: 4.51" (114.55mm)
Outer Diameter of Tube 1.16" (29.46mm)
Inner Diameter of Tube 1.08" (27.43mm)

Packaging: Racked (cardboard fiber racks)
25 tips per rack (with attached caps)
2 racks per pack (50 tubes)
10 packs per case (500 tubes total)
500 tubes/caps

Autoclavable. Freezable. Tubes can be centrifuged to 20,000 RCF.

MicroBiological testing:

Certified testing is performed lot by lot. These tests include Nuclease (RNase and DNase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Sterility:

Certified as sterile with a sterility assurance level (SAL) of 10^{-6}