



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

## SuperClear® 12x75 mm Culture Tubes with Plug Caps, Polystyrene, 25 per Bag, Sterile

Fluid Measure	mL	Freezable	No
Fluid Volume	5	Freeze Temp	-80°C
Height (mm)	75	Autoclavable	No
Width (mm)	12	Sterile	Yes
Length (mm)		Sterile Processing	Certified as sterile with a sterility assurance level (SAL) of 10 <sup>-6</sup>
Pack Qty	25	Centrifugation	Yes
Case Qty	500	Centrifugation Speed (RCF)	1500
Package Type	Bag	Shelf Life (Months)	36
Material	Polystyrene	Certification	
Cap Material	Polyethylene	Regulatory	

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing: ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

