



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

## SuperClear® 13x100 mm Culture Tubes, Polypropylene, 125 per Bag

Fluid Measure	mL	Freezable	Yes
Fluid Volume	8	Freeze Temp	-80°C
Height (mm)	100	Autoclavable	121°C 15min
Width (mm)	13	Sterile	No
Length (mm)		Sterile Processing	NA
Pack Qty	125	Centrifugation	Yes
Case Qty	1000	Centrifugation Speed (RCF)	3000
Package Type	Bag	Shelf Life (Months)	48
Material	Polypropylene	Certification	
Cap Material		Regulatory	

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing: ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

