



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

SuperClear® 13x100 mm Culture Tubes, Polypropylene, 125 per Bag, Sterile

Fluid Measure	mL	Freezable	No
Fluid Volume	8	Freeze Temp	-80°C
Height (mm)	100	Autoclavable	121°C 15min
Width (mm)	13	Sterile	Yes
Length (mm)		Sterile Processing	Certified as sterile with a sterility assurance level (SAL) of 10 ⁻⁶
Pack Qty	125	Centrifugation	Yes
Case Qty	1000	Centrifugation Speed (RCF)	3000
Package Type	Bag	Shelf Life (Months)	36
Material	Polypropylene	Certification	
Cap Material		Regulatory	

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing: ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

