



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

Dual Position Caps for 17x100 mm Culture Tubes, in Bags

Fluid Measure		Freezable	Yes
Fluid Volume		Freeze Temp	-80°C
Height (mm)	18.79	Autoclavable	No
Width (mm)	21.84	Sterile	No
Length (mm)		Sterile Processing	NA
Pack Qty	1000	Centrifugation	NA
Case Qty	1000	Centrifugation Speed (RCF)	
Package Type	Bag	Shelf Life (Months)	48
Material	Polyethylene	Certification	
Cap Material		Regulatory	

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing: Endotoxin Free (non-pyrogenic)

