



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

SuperClear® 17x100 mm Culture Tubes with Separate Dual Position Caps, Polypropylene, in Bags

Fluid Measure	mL	Freezable	Yes
Fluid Volume	14	Freeze Temp	-80°C
Height (mm)	100	Autoclavable	121°C 15min
Width (mm)	17	Sterile	No
Length (mm)		Sterile Processing	NA
Pack Qty	1000	Centrifugation	Yes
Case Qty	1000	Centrifugation Speed (RCF)	3000
Package Type	Bag	Shelf Life (Months)	48
Material	Polypropylene	Certification	
Cap Material	Polyethylene	Regulatory	

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing: ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

