



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

## SuperClear® 17x100 mm Culture Tubes with Dual Position Caps, Polypropylene, 25 per Bag, Sterile

Fluid Measure	mL	Freezable	Yes
Fluid Volume	14	Freeze Temp	-80°C
Height (mm)	100	Autoclavable	121°C 15min
Width (mm)	17	Sterile	Yes
Length (mm)		Sterile Processing	Certified as sterile with a sterility assurance level (SAL) of 10 <sup>-6</sup>
Pack Qty	25	Centrifugation	Yes
Case Qty	500	Centrifugation Speed (RCF)	3000
Package Type	Bag	Shelf Life (Months)	36
Material	Polypropylene	Certification	
Cap Material	Polyethylene	Regulatory	

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing: ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

