

## Product Description

**Catalog Number:** 3631-870-000-9

**Product Description:** 0.5ml Micro Tube with Screw Cap  
Polypropylene. Clear. Non-Sterile

**Product Dimensions:**

Length of Tube: 1.57" (39.87mm)  
Outer Diameter: 0.41" (10.41mm)  
Inner Diameter: 0.33" (8.38mm)

**Packaging:** Bagged  
500 tubes per bag  
10 bags per case (5000 tubes total)

\*Screw caps are in separate bags

Freezable. Tubes can be centrifuged to 20,000 RCF.

Autoclavable. (Remember to loosen the caps of tubes to prevent collapse when autoclaving)

**MicroBiological testing:**

Certified testing is performed lot by lot. These tests include Nuclease (RNase and DNase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

**Performance Testing:**

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.