



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

## SuperClear® 0.5 mL Freestanding Screw Cap Microcentrifuge Tubes with Elastomeric Caps, Sterile

Fluid Measure	mL	Freezable	Yes
Fluid Volume	0.5	Freeze Temp	-80°C
Height (mm)	45.45	Autoclavable	121°C 15min
Width (mm)	10.4	Sterile	Yes
Length (mm)		Sterile Processing	Certified as sterile with a sterility assurance level (SAL) of 10 <sup>-6</sup>
Pack Qty	500	Centrifugation	Yes
Case Qty	5000	Centrifugation Speed (RCF)	20000
Package Type	Resealable Bags	Shelf Life (Months)	36
Material	Polypropylene	Certification	
Cap Material	Polypropylene	Regulatory	

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing: ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

