

Dimensions are inches. Do not scale




Image is not actual size. Colors may differ.

PurePlus® 0.1 mL 96 Well Low Profile PCR Plates with Full Skirt for Popular Thermocyclers, Natural Color

Fluid Measure	mL	Freezable	Yes
Fluid Volume	0.1	Freeze Temp	-80°C
Height (mm)	15.5	Autoclavable	No
Width (mm)	85.5	Sterile	No
Length (mm)	127.8	Sterile Processing	NA
Pack Qty	10	Centrifugation	NA
Case Qty	100	Centrifugation Speed (RCF)	
Package Type	Box	Shelf Life (Months)	48
Material	Polypropylene	Certification	
Cap Material		Regulatory	

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing: ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

This Technical Data Sheet is provided as a service and for documentation only. Labcon North America retains all rights to the information, procedures, and results of the material contained hereon as proprietary information. Any copying and distribution of this technical data without prior written authorization is prohibited by law. The use of this information and documentation does not discharge the user from their responsibility to test and use this product appropriately. This item is for use solely by laboratory professionals. This document is automatically generated and is © Copyright 04-Mar-2022 . All Rights Reserved. Technical Data Sheet Version 1/2020

DOMINIQUE DUTSCHER SAS