


Technical Product Data

Technical Data Sheet

Labcon Product Number 3996-520-000-9



Dimensions are millimeters. Do not scale

Image is not actual size. Colors may differ.

Reusable Silicone Sealing Mat for Square Hole Plates, Autoclavable

Fluid Measure		Freezable	Yes
Fluid Volume		Freeze Temp	-80°C
Height (mm)	4.53	Autoclavable	121°C 15min
Width (mm)	76.7	Sterile	
Length (mm)	119.2	Sterile Processing	Non-sterile
Pack Qty	5	Centrifugation	NA
Case Qty	50	Centrifugation Speed (RCF)	
Package Type	Box	Shelf Life (Months)	48
Material	Silicone	Certification	
Cap Material		Regulatory	

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing: Endotoxin Free (non-pyrogenic)



This Technical Data Sheet is provided as a service and for documentation only. Labcon North America retains all rights to the information, procedures, and results of the material contained herein as proprietary information. Any copying and distribution of this technical data without prior written authorization is prohibited by law. The use of this information and documentation does not discharge the user from their responsibility to test and use this product appropriately. This item is for use solely by laboratory professionals. This document is automatically generated and is © Copyright 25-Jan-2021 . All Rights Reserved. Technical Data Sheet Version 1/2020