



Solution Bottles



- Compatible with bottle top filters or for use as solution storage containers (Bottle Top Filters sold separately)

- Solution bottle neck size is 45mm
- Curved bottle design for easy handling




SOLUTION BOTTLES			S, R/D F, NP, USP, PS	
Part No.	Product Description	Working volume	Pack Description	Qty/Case
229781	150mL Solution Bottle	150mL	1/Bag	24
229782	250mL Solution Bottle	250mL	1/Bag	24
229784	500mL Solution Bottle	500mL	1/Bag	24
229785	1000mL Solution Bottle	1000mL	1/Bag	24

500mL Square Media Bottle

- Ideal for media prep and storage of laboratory reagents and biological buffers
- Square design saves valuable shelf space and is ideal for stacking, such as in a refrigerator
- Unique individually bagged packaging allows opening one at a time
- PETG material reduces gas permeability

- Molded graduations on two sides
- Flat surface polyethylene (HDPE) caps are ideal for marking



SQUARE MEDIA BOTTLE			S, R/D F, NP, USP, PETG	
Part No.	Product Description	Working volume	Pack Description	Qty/Case
229788	500mL Media Bottle, Square	500mL	1/Bag	24