

## Solution Bottles



- Compatible with bottle top filters or for use as solution storage containers (Bottle Top Filters sold separately)

- 1 Solution bottle neck size is 45mm
- 2 Curved bottle design for easy handling



### SOLUTION BOTTLES

S, R/D F, NP, USP, PS

Part No.	Product Description	Working volume	Pack Description	Qty/Case
229781	150mL Solution Bottle	150mL	1/Bag	24
229782	250mL Solution Bottle	250mL	1/Bag	24
229784	500mL Solution Bottle	500mL	1/Bag	24
229785	1000mL Solution Bottle	1000mL	1/Bag	24

## 500mL Square Media Bottle

- Ideal for media prep and storage of laboratory reagents and biological buffers
- Square design saves valuable shelf space and is ideal for stacking, such as in a refrigerator
- Unique individually bagged packaging allows opening one at a time
- PETG material reduces gas permeability

- 1 Molded graduations on two sides
- 2 Flat surface polyethylene (HDPE) caps are ideal for marking



### SQUARE MEDIA BOTTLE

S, R/D F, NP, USP, PETG

Part No.	Product Description	Working volume	Pack Description	Qty/Case
229788	500mL Media Bottle, Square	500mL	1/Bag	24