

Disposable Aspirating Pipet Specification

Product: 1mL Aspirating Pipet

Features:

- Easy view, non-printed pipet
- Same dimensions and quality of serological pipet.
- Drawn tip
- No plug
- Pipets are available individually wrapped or in bulk pack
- Pipets are supplied sterile, by gamma irradiation

Material: Clear virgin polystyrene, USP Class VI, Non-Pyrogenic

Specifications: A - Overall Length – 270mm (10.6in) +/-1.0mm
B – Tip ID – 1.6mm (0.06in) +/-0.2mm
C – Mouth Piece OD – 4.5mm (0.18in) +/-0.15mm
D – Tube OD - 4.5mm (0.18in) +/-0.15mm

Volume: 1.5mL total volume

Product Offering: Celltreat P/N 229261 - Individually Paper/Poly Wrapped, 100/Bag, 400/Case
Celltreat P/N 229251 - Bulk Packed, 50/Bag, 400/Case

