

FICHE TECHNIQUE / DATASHEET

Référence	978361 978362
Désignation	Flacons de culture / Culture flasks: 12,5 cm ²



Description :

Reference / Part Number	Description
978361	12.5cm ² Tissue Culture Flask - Plug Seal Cap, Sterile
978362	12.5cm ² Tissue Culture Flask - Vent Cap, Sterile

- ✓ Uniform hydrophilic surface
Surface hydrophile uniforme
- ✓ Vent caps feature 0.22µm hydrophobic membrane **(978362 only)**
Capuchon ventilé à membrane hydrophobe 0,22µm
- ✓ Canted neck design for easy access
Col incliné pour un accès facile
- ✓ Packaged sterile (gamma irradiated) in zip-closure bags
Emballés stériles (irradiés gamma) en sachets zippés
- ✓ Tissue culture treated surface
Surface traitée pour la culture cellulaire
- ✓ Capacity : 25 ml
- ✓ Material : Polystyrene

Unité de vente (Sales unit) :

200 Flasks